

WIDEX RC-DEX

REMOTE CONTROL



MAIN FEATURES

- Program shift
- Volume control
- LED indicating transmission
- Battery replaceable by user
- Battery life approx. 12 months
- Attachable to key ring
- Available in silver or black color

RADIO SPECIFICATIONS

Frequency (range)	10.6 MHz (10.2-11.0 MHz)
Bandwidth (-15 dB)	660 kHz
Channel	Single channel radio
Modulation	FSK
Rated output power	21 nW/-46.8 dBm
Magnetic field strength	-13 dBμA/m @ 10 m
Duty cycle	0.1 %
	Simplex
Conforms to the following standards/approvals	• EN/IEC 60950-1 • EN 50371 • EN 300 330-2 • EN 301 489-1 • EN 301 489-3 • FCC Part 15 C • FCC ID: TTY-RCDEX • CE Mark

PHYSICAL PROPERTIES

Dimensions:	2.0/1.2/0.6 in (h/w/d)
Weight:	0.60 oz
Battery type:	Panasonic CR2032 lithium, 3 Volts/200 mAh


ACCESSORIES

Battery
Key ring/string

RC-DEX

OVERVIEW

1. LED
2. Volume up
3. Program toggle
4. Volume down



1. Battery cover
2. Eye for key ring/string attachment

