


WIDEX UNI-DEX

AUDIO CONTROL



MAIN FEATURES

- Hands-free listening solutions for external audio devices
- Connects to mobile phones, computers, tablets
- Connects to all devices with a 3.5 mm jack output
- Use hearing aids as headset
- Room off function
- Automatic start of streaming
- Long battery life, 40 hours streaming/ 4 months standby
- Simple and easy to use
- Fast battery charging, 1 hour
- Wireless audio streaming in high quality to hearing aids with WidexLink

RADIO SPECIFICATIONS

The UNI-DEX contains a radio transmitter	
Radio transmitter parameters:	
Frequency (range)	10.6 MHz (10.2-11.0 MHz)
Bandwidth (-15 dB)	660 kHz
Channel	Single channel radio
Modulation	FSK
Rated output power	3.02 pW / -85.2 dBm (ERP effective radiated power)
Magnetic field strength	-26.9 dBμA/m @ 10 m
Duty cycle	100%
	Simplex capability - transmitting only
Conforms to the following standards/approvals	EN/IEC 60950-1 EN 62479 EN 300 330-2 EN 301 489-1 EN 301 489-3 EN/IEC 60601-1-2 FCC Part 15 C IC RSS-210 CE mark


PHYSICAL PROPERTIES

Dimensions:	1.65/0.79/0.35 inch
Weight:	0.61 oz
Battery type:	Lithium-ion Polymer
Audio cable:	29.53 inch
Neck-loop:	3.94 inch & 25.98 inch plus connectors - 29.92 inch closed

UNI-DEX

OVERVIEW

1. Push button/answer button
2. Room off button
3. Neck loop
4. Microphone
5. Charging indicator
6. Charger connector
7. Jack plug



5. Charging indicator
6. Charger connector

