

WIDEX COM-DEX

WIRELESS STREAMER



MAIN FEATURES

- Hands-free voice call for most Bluetooth enabled phones
- Music streaming for most Bluetooth enabled phones
- Smartphone App for hearing aid remote control
- Hearing aids can be used as a headset/hands-free solution
- High-fidelity audio streaming to hearing aids with WidexLink
- Room Off function
- Automatic streaming
- Constant streaming up to 8 hours/8 days standby*
- Battery charging 3 hours
- Stylish colors
- Simple and easy to use

RADIO SPECIFICATIONS

Radio transmitter parameters:

	Bluetooth	WidexLink™
Version	2.1 + EDR	N/A
Profiles	A2DP, AVRCP, HFP/HSP	N/A
Frequency (range)	2.4-2.5 GHz	10.6 MHz (10.2-11.0 MHz)
Bandwidth	1 MHz	660 kHz (-15 dB)
Channel	5 logical channels	Single channel radio
Modulation	FHSS/GFSK, $\pi/4$ DPSK, 8DPSK	FSK
Radiated output power	2.5 mW ($\approx +4.0$ dBm)	3.02 pW / -85.2 dBm (ERP effective radiated power)
Magnetic field strength	N/A	-26.9 dB μ A/m @ 11 yd
Range	Up to 11 yd to mobile phone	Up to 1 ft from neck loop to hearing aid(s)
Duty cycle	Max 83 % Time division duplex (TDD)	100% Simplex or semi duplex
Conforms to the following standards/approvals	EN/IEC 60950-1 EN 62479 EN 300 328 EN 300 330-2 EN 301 489-3 EN 301 489-17 EN/IEC 60601-1-2 FCC Part 15 C IC RSS-210 IC RSS-247 CE mark	

PHYSICAL PROPERTIES

Dimensions:	2.76/0.79/0.31 inch (l/w/h)
Neck loop:	30 inch
Weight:	0.95 oz
Battery type:	Lithium-ion Polymer, rechargeable

*12 hours per day

COM-DEX

OVERVIEW

1. Push button/Answer button
 2. Room Off button
 3. Neck loop
 4. Microphone
 5. LED Indicator
 6. Charger connector/MICRO USB
 7. On/off switch by opening/closing the neck loop
-

