

WIDEX CALL-DEX EASY PLUG AND PLAY



MAIN FEATURES

- Discreet and compact design
- Direct streaming via WidexLink™
- Superior sound quality
- Plug and play in most phones (3.5 mm jack)
- Up to 80 hours of constant streaming
- Fully automatic (no buttons)

RADIO SPECIFICATIONS

The CALL-DEX contains a radio transmitter	
Radio transmitter parameters:	WidexLink™
Frequency (range)	10.6 MHz (10.2-11.0 MHz)
Bandwidth (-15 dB)	660 kHz (-15 dB)
Channel	Single channel radio
Modulation	FSK
Radiated output power	7 pW (≈ -81.0 dBm)
Magnetic field strength	-45 dBμA/m @ 10 m
Duty cycle	100%
Conforms to the following standards/approvals	EN/IEC 60950-1 EN 62479 EN 300 330-2 EN 301 489-1 EN 301 489-3 EN/IEC 60601-1-2 FCC Part 15 C IC RSS-210 CE mark

PHYSICAL PROPERTIES

Dimensions:	0.87/0.35/0.47 inch. (l/w/h)
Weight:	0.31 oz
Battery type:	ZA10

CALL-DEX

OVERVIEW

1. Jack output
 2. Battery lid
-

